

Kallikrein

Chemical Properties

CAS No. : 9001-01-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Kallikrein is capable of forming the kallikrenase kalinin system (KKS) in plasma and tissues, producing bradykinin and kalin peptides , respectively [1] . |
| Targets(IC50) | Others,Endogenous Metabolite |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481