

Peginterferon beta-1a

Chemical Properties

CAS No. : 1211327-92-2

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Peginterferon beta-1a (Peginterferon β -1a), the inaugural pegylated interferon beta-1a derivative, promotes apoptosis in cancer cells and exhibits anti-tumor properties in nude mice models. It is utilized in the study of cancer and relapsing multiple sclerosis (RMS) [1].
Targets(IC50)	Apoptosis,Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481