

3-Hydroxysebacic acid

Chemical Properties


CAS No. : 73141-46-5

Formula: C₁₀H₁₈O₅

Molecular Weight: 218.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Hydroxysebacic acid, an endogenous metabolite found in urine, is utilized in the study of Medium Chain Acyl-CoA Dehydrogenase Deficiency (MCADD) [1] [2].
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5819 mL	22.9095 mL	45.819 mL
5 mM	0.9164 mL	4.5819 mL	9.1638 mL
10 mM	0.4582 mL	2.291 mL	4.5819 mL
50 mM	0.0916 mL	0.4582 mL	0.9164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481