


## Normetanephine

## Chemical Properties

CAS No. :	97-31-4
Formula:	C <sub>9</sub> H <sub>13</sub> NO <sub>3</sub>
Molecular Weight:	183.2044
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Normetanephine, an endogenous metabolite found in cerebrospinal fluid, has applications in hypertension research [1] [2].
Targets(IC50)	Others,Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4585 mL	27.2926 mL	54.5852 mL
5 mM	1.0917 mL	5.4585 mL	10.917 mL
10 mM	0.5459 mL	2.7293 mL	5.4585 mL
50 mM	0.1092 mL	0.5459 mL	1.0917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481