

## PROTAC BTK Degradar-2

## Chemical Properties

CAS No. : 2250382-66-0


Formula: C47H54F2N8O13

Molecular Weight: 976.97

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTAC BTK Degradar-2, a potent degrader of BTK through the PROTAC mechanism, effectively diminishes BTK protein levels [1].
Targets(IC50)	Others,BTK,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0236 mL	5.1179 mL	10.2357 mL
5 mM	0.2047 mL	1.0236 mL	2.0471 mL
10 mM	0.1024 mL	0.5118 mL	1.0236 mL
50 mM	0.0205 mL	0.1024 mL	0.2047 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481