

DSPE-mPEG-2000

Chemical Properties

CAS No. : 156543-00-9

Formula: C₄₅H₈₈N₀O₁₁P

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | DSPE-mPEG-2000, a PEG-based phospholipid, facilitates the synthesis of liposomes for cancer agent delivery [1] [2]. |
| Targets(IC50) | Others,Liposome |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481