

DSPE-PEG-PDP

Chemical Properties

CAS No. : 474922-24-2
Formula: (C₂H₄O)_nC₅₂H₉₄N₃O₁₁P₅.H₃N
Molecular Weight:
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-PDP, a phospholipid-polyethylene glycol (PEG) conjugate, is utilized in drug delivery applications [1].
Targets(IC50)	Others,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481