

## DSPE-PEG-Carboxylic Acid

### Chemical Properties

CAS No. : 474922-20-8

Formula: C<sub>46</sub>H<sub>91</sub>N<sub>2</sub>O<sub>13</sub>P

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DSPE-PEG-Carboxylic Acid, a phospholipid polyethylene glycol (PEG) conjugate, facilitates the targeted delivery of agents and genes [1].
Targets(IC50)	Others,Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481