

## C8 PEG-Ceramide

### Chemical Properties

CAS No. : 212116-76-2  
Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>31</sub>H<sub>57</sub>N<sub>06</sub>  
Molecular Weight:  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	C8 PEG-Ceramide, a lipid product, possesses the capability to synthesize lipid bilayer carriers, making it applicable for drug delivery [1] [2].
Targets(IC50)	Others,Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481