

DSPE-PEG-Cyanur

Chemical Properties

CAS No. : 1246304-74-4

Formula: C₄₉H₉₃Cl₂N₆O₁₁P

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-Cyanur, a PEG derivative with a cyanur (cyanuric) functional group, is utilized for the PEGylation of proteins under mild basic conditions and serves as a component for nanostructured lipid carriers [1] [2].
Targets(IC50)	Others,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481