

mPEG-DSPE ammonium

## Chemical Properties

CAS No. : 474922-77-5

Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>43</sub>H<sub>84</sub>N<sub>10</sub>P.H<sub>3</sub>N

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	18:0 PEG-PE, a polyethyleneglycol/phosphatidyl-ethanolamine (PEG-PE) conjugate, serves in drug delivery applications [1].
Targets(IC50)	Others,Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481