


DNA crosslinker 1 dihydrochloride

Chemical Properties

CAS No. :	2761734-27-2
Formula:	C15H22Cl2N8O
Molecular Weight:	401.29
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	DNA Crosslinker 1 (dihydrochloride) serves as a potent binder for the DNA minor groove, exhibiting a DNA binding affinity (ΔT_m) of 1.1°C. This compound is utilized in research focused on anticancer therapies [1].
Targets(IC50)	Others,DNA Alkylator/Crosslinker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.492 mL	12.4598 mL	24.9196 mL
5 mM	0.4984 mL	2.492 mL	4.9839 mL
10 mM	0.2492 mL	1.246 mL	2.492 mL
50 mM	0.0498 mL	0.2492 mL	0.4984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481