

## PROTAC SOS1 degrader-2

## Chemical Properties

CAS No. : 2913176-81-3


Formula: C57H76ClFN10O4S

Molecular Weight: 1051.79

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTAC SOS1 degrader-2, a potent degrader of PROTAC SOS1, effectively reduces the expression of pERK and RAS-GTP in a dose-dependent manner. It significantly inhibits tumor growth in vivo [1].
Targets(IC50)	Others,PROTACs,Ras

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9508 mL	4.7538 mL	9.5076 mL
5 mM	0.1902 mL	0.9508 mL	1.9015 mL
10 mM	0.0951 mL	0.4754 mL	0.9508 mL
50 mM	0.019 mL	0.0951 mL	0.1902 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481