

CDK4/6-IN-11

Chemical Properties

CAS No. : 2139329-47-6


Formula: C43H49N11O7

Molecular Weight: 831.92

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CDK4/6-IN-11 is a potent PROTAC CDK4/6 degrader.
Targets(IC50)	Others,CDK,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.202 mL	6.0102 mL	12.0204 mL
5 mM	0.2404 mL	1.202 mL	2.4041 mL
10 mM	0.1202 mL	0.601 mL	1.202 mL
50 mM	0.024 mL	0.1202 mL	0.2404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481