

Sozinibercept

Chemical Properties

CAS No. : 2568358-31-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sozinibercept (OPT 302; VGX-300) is a soluble form of Vascular Endothelial Growth Factor Receptor-3 (VEGFR-3) that potently inhibits the activity of Vascular Endothelial Growth Factor-C and -D (VEGF-C/D), key proangiogenic factors involved in promoting angiogenesis and vascular leakage. Additionally, Sozinibercept has been shown to inhibit diabetic retinal edema in rats [1] [2] [3].
Targets(IC50)	VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481