

Ogalvibart

Chemical Properties

CAS No. : 2599039-60-6

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ogalvibart (C-135-LS), a humanized anti-SARS-CoV-2 IgG1 type monoclonal antibody, specifically targets and binds to the spike (S) glycoprotein receptor-binding domain (RBD) of the SARS-CoV-2 virus. When combined with C144LS in a 1:1 ratio, Ogalvibart demonstrates significant preventive activity and effectively prevents the progression of COVID-19 in a rhesus monkey disease model [1].
Targets(IC50)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481