

m7Gpppm6AmpG ammonium

Chemical Properties

CAS No. :

Formula: C33H45N15O24P4.xNH3

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	M7Gpppm6AmpG ammonium, a trinucleotide mRNA 5' cap analog, facilitates in vitro RNA synthesis [1].
Targets(IC50)	Others,DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481