


## PROTAC VEGFR-2 degrader-2

## Chemical Properties

CAS No. :	2353417-85-1
Formula:	C50H56FN9O6S
Molecular Weight:	930.1
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	PROTAC VEGFR-2 degrader-2 (PROTAC-4), a specific degrader of VEGFR-2, demonstrates minimal inhibition of VEGFR-2 (IC <sub>50</sub> > 1 μM) and low anti-proliferative activity towards EA.hy926 cells (IC <sub>50</sub> > 100 μM) [1].
Targets(IC <sub>50</sub> )	Others,PROTACs,VEGFR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0752 mL	5.3758 mL	10.7515 mL
5 mM	0.215 mL	1.0752 mL	2.1503 mL
10 mM	0.1075 mL	0.5376 mL	1.0752 mL
50 mM	0.0215 mL	0.1075 mL	0.215 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481