


STAT3 degrader-1

Chemical Properties

CAS No. :	2497585-16-5
Formula:	C ₅₈ H ₆₃ F ₅ N ₉ O ₁₂ PS
Molecular Weight:	1236.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	STAT3 Degrader-1 (Compound 295) is a potent degrader of STAT3, utilized in anticancer research [1].
Targets(IC50)	Others,STAT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8089 mL	4.0447 mL	8.0893 mL
5 mM	0.1618 mL	0.8089 mL	1.6179 mL
10 mM	0.0809 mL	0.4045 mL	0.8089 mL
50 mM	0.0162 mL	0.0809 mL	0.1618 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481