

Non-PEGylated/naked ARC1905

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ARC1905, a non-PEGylated/naked anti-C5 RNA aptamer, functions by inhibiting the cleavage of complement factor 5 (C5) into its components, C5a and C5b, without PEG conjugation. This compound is currently being researched for its potential application in the treatment of age-related macular degeneration (AMD).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481