

EDA-DA

Chemical Properties


CAS No. : 87156-00-1

Formula: C₈H₁₂N₂O₃

Molecular Weight: 184.19

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	EDA-DA, an N-terminally tagged dipeptide probe, serves as a labeling agent for bacterial Peptidoglycan (PG), a fundamental component of the cell walls in most bacteria crucial for cellular division, shape retention, and hydrostatic pressure maintenance [1].
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4292 mL	27.1459 mL	54.2918 mL
5 mM	1.0858 mL	5.4292 mL	10.8584 mL
10 mM	0.5429 mL	2.7146 mL	5.4292 mL
50 mM	0.1086 mL	0.5429 mL	1.0858 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481