

Elebsiran

Chemical Properties

CAS No. :	2648009-64-5
Formula:	C485H658F9N158O294P39S6
Molecular Weight:	14976.62
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

Biological Description

Description	Elebsiran is an antiviral agent [1] .
Targets(IC50)	Others,HBV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0668 mL	0.3339 mL	0.6677 mL
5 mM	0.0134 mL	0.0668 mL	0.1335 mL
10 mM	0.0067 mL	0.0334 mL	0.0668 mL
50 mM	0.0013 mL	0.0067 mL	0.0134 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481