


Ilaprazole sulfone

Chemical Properties

CAS No. :	172152-37-3
Formula:	C19H18N4O3S
Molecular Weight:	382.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ilaprazole sulfone, primarily metabolized by CYP3A4/5, constitutes the chief metabolite of Ilaprazole (IY-81149), an orally active proton pump inhibitor [1].
Targets(IC50)	Others,Drug Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6148 mL	13.0739 mL	26.1479 mL
5 mM	0.523 mL	2.6148 mL	5.2296 mL
10 mM	0.2615 mL	1.3074 mL	2.6148 mL
50 mM	0.0523 mL	0.2615 mL	0.523 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481