


Abiraterone N-oxide

Chemical Properties

CAS No. : 2378463-76-2
 Formula: C₂₄H₃₁NO₂
 Molecular Weight: 365.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Abiraterone N-oxide, a metabolite of Abiraterone, serves as a potent and irreversible CYP17A1 inhibitor with antiandrogen activity [1].
Targets(IC50)	Others,Cytochromes P450,Drug Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7359 mL	13.6795 mL	27.359 mL
5 mM	0.5472 mL	2.7359 mL	5.4718 mL
10 mM	0.2736 mL	1.368 mL	2.7359 mL
50 mM	0.0547 mL	0.2736 mL	0.5472 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481