

TRPV1 activator-2

Chemical Properties


CAS No. : 17514-11-3

Formula: C₁₉H₂₉NO₃

Molecular Weight: 319.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TRPV1 activator-2 (Compound 9), a capsaicin head analog, engages in targeted interactions with channel residues at the lipid-water interface [1].
Targets(IC50)	Others,TRP/TRPV Channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1305 mL	15.6524 mL	31.3048 mL
5 mM	0.6261 mL	3.1305 mL	6.261 mL
10 mM	0.313 mL	1.5652 mL	3.1305 mL
50 mM	0.0626 mL	0.313 mL	0.6261 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481