


## ER degrader 5

## Chemical Properties

CAS No. :	2913192-47-7
Formula:	C <sub>26</sub> H <sub>18</sub> F <sub>2</sub> O <sub>4</sub> S
Molecular Weight:	464.48
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ER Degrader 5 is a potent compound that degrades the estrogen receptor (ER), exhibits anti-proliferation activity, and has potential applications in breast cancer research [1].
Targets(IC50)	Estrogen Receptor/ERR,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1529 mL	10.7647 mL	21.5295 mL
5 mM	0.4306 mL	2.1529 mL	4.3059 mL
10 mM	0.2153 mL	1.0765 mL	2.1529 mL
50 mM	0.0431 mL	0.2153 mL	0.4306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481