


## TAMRA alkyne, 6-isomer

### Chemical Properties

CAS No. :	1352649-44-5
Formula:	C <sub>28</sub> H <sub>25</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	467.52
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



### Biological Description

Description	TAMRA alkyne, 6-isomer, a variant of the xanthene dye TAMRA, known for its orange emission, serves as a linker specifically designed for oligonucleotide labeling and amino acid sequencing. This isomer incorporates alkyne groups, enabling copper-catalyzed Click Chemistry reactions with azides [1].
Targets(IC50)	Others, Autophagy

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1389 mL	10.6947 mL	21.3895 mL
5 mM	0.4278 mL	2.1389 mL	4.2779 mL
10 mM	0.2139 mL	1.0695 mL	2.1389 mL
50 mM	0.0428 mL	0.2139 mL	0.4278 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481