


Pirnabine

Chemical Properties

CAS No. :	19825-63-9
Formula:	C ₁₉ H ₂₄ O ₃
Molecular Weight:	300.39
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Pirnabine, a cannabinoid receptor ligand, is utilized in glaucoma research [1].
Targets(IC50)	Cannabinoid Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.329 mL	16.645 mL	33.2901 mL
5 mM	0.6658 mL	3.329 mL	6.658 mL
10 mM	0.3329 mL	1.6645 mL	3.329 mL
50 mM	0.0666 mL	0.3329 mL	0.6658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481