

2-(n-Propylidene hydrazino) adenosine

Chemical Properties

CAS No. :

Formula: C₁₃H₁₉N₇O₄

Molecular Weight: 337.33

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-(n-Propylidene hydrazino) adenosine is a purine nucleoside analog. Purine nucleoside analogs have broad antitumor activity targeting indolent lymphoid malignancies. Anticancer mechanisms in this process rely on inhibition of DNA synthesis, induction of apoptosis , etc [1] .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9645 mL	14.8223 mL	29.6446 mL
5 mM	0.5929 mL	2.9645 mL	5.9289 mL
10 mM	0.2964 mL	1.4822 mL	2.9645 mL
50 mM	0.0593 mL	0.2964 mL	0.5929 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481