

## Bamosiran

## Chemical Properties

CAS No. : 1337968-84-9

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Bamosiran, a small interfering RNA (siRNA), specifically targets the $\beta$ -adrenergic receptor 2 to reduce intraocular pressure.
Targets(IC50)	Others, Adrenergic Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481