


RDR 02308

Chemical Properties

CAS No. : 4155-82-2
 Formula: C₁₉H₁₅N₃O₃
 Molecular Weight: 333.34
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	RDR 02308 is a TLR4-MyD88 binding inhibitor that inhibits full-length β -lactamase [1].
Targets(IC50)	Others,Antibacterial,TLR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9999 mL	14.9997 mL	29.9994 mL
5 mM	0.600 mL	2.9999 mL	5.9999 mL
10 mM	0.300 mL	1.500 mL	2.9999 mL
50 mM	0.060 mL	0.300 mL	0.600 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481