


UK-2A

Chemical Properties

CAS No. : 167173-85-5
 Formula: C₂₆H₃₀N₂O₉
 Molecular Weight: 514.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	UK-2A (Antibiotic UK 2A), a potent antifungal antibiotic, demonstrates significant antifungal activity [1] [2].
Targets(IC50)	Others,Antibiotic,Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9436 mL	9.7178 mL	19.4356 mL
5 mM	0.3887 mL	1.9436 mL	3.8871 mL
10 mM	0.1944 mL	0.9718 mL	1.9436 mL
50 mM	0.0389 mL	0.1944 mL	0.3887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481