

alpha-2-3,6-sialidase (BiNanH2)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alpha-2-3,6-sialidase (BiNanH2), a sialyltransferase commonly utilized in biochemical research, catalyzes the transfer of $\alpha$ -2,3/2,6-sialyl from CMP-Neu5Ac to galactoside acceptors [1].
Targets(IC50)	Others,Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481