

alpha-2,3-Sialyltransferase (PmST3)

Chemical Properties

CAS No. : 71124-51-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Alpha-2,3-Sialyltransferase (PmST3) is a beta-galactoside that catalyzes the transfer of sialic acid to carbohydrate groups of glycoproteins and glycolipids [1].
Targets(IC50)	Others,Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481