


Orotic acid potassium

Chemical Properties

CAS No. :	24598-73-0
Formula:	C ₅ H ₃ KN ₂ O ₄
Molecular Weight:	194.19
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Orotic acid (potassium salt) serves as a biochemical reagent, applicable as both a biological material and an organic compound in life science research.
Targets(IC50)	Others,Endogenous Metabolite,DNA/RNA Synthesis

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1496 mL	25.748 mL	51.496 mL
5 mM	1.0299 mL	5.1496 mL	10.2992 mL
10 mM	0.515 mL	2.5748 mL	5.1496 mL
50 mM	0.103 mL	0.515 mL	1.0299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481