


Monaschromone

Chemical Properties

CAS No. :	1338576-70-7
Formula:	C ₁₁ H ₁₂ O ₄
Molecular Weight:	208.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Monaschromone, a polyketide metabolite, effectively suppresses the growth of <i>B. cinerea</i> , <i>A. solani</i> , <i>M. oryzae</i> , and <i>G. saubinettii</i> , demonstrating minimum inhibitory concentration (MIC) values between 6.25 and 12.5 μ M [1].
Targets(IC50)	Others,Endogenous Metabolite,Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8028 mL	24.0142 mL	48.0284 mL
5 mM	0.9606 mL	4.8028 mL	9.6057 mL
10 mM	0.4803 mL	2.4014 mL	4.8028 mL
50 mM	0.0961 mL	0.4803 mL	0.9606 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481