


## Bipolamine G

## Chemical Properties

CAS No. :	1357824-71-5
Formula:	C <sub>21</sub> H <sub>28</sub> N <sub>2</sub> O <sub>4</sub>
Molecular Weight:	372.46
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bipolamine G is an antibacterial polyketide alkaloid [1].
Targets(IC50)	Others,Endogenous Metabolite,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6849 mL	13.4243 mL	26.8485 mL
5 mM	0.537 mL	2.6849 mL	5.3697 mL
10 mM	0.2685 mL	1.3424 mL	2.6849 mL
50 mM	0.0537 mL	0.2685 mL	0.537 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481