


Blestrin D

Chemical Properties

CAS No. :	145701-19-5
Formula:	C ₃₀ H ₂₄ O ₆
Molecular Weight:	480.51
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Blestrin D, a potent mixed-type butyrylcholinesterase (BChE) inhibitor with an IC ₅₀ value of 8.1 μM, is isolated from <i>Bletilla striata</i> and has applications in Alzheimer's disease (AD) research [1].
Targets(IC50)	Others,Cholinesterase (ChE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0811 mL	10.4056 mL	20.8112 mL
5 mM	0.4162 mL	2.0811 mL	4.1622 mL
10 mM	0.2081 mL	1.0406 mL	2.0811 mL
50 mM	0.0416 mL	0.2081 mL	0.4162 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481