


Justicidin C

Chemical Properties

CAS No. :	17803-12-2
Formula:	C ₂₂ H ₁₈ O ₇
Molecular Weight:	394.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Justicidin C, a lignan compound, exhibits potent antiviral properties, particularly against vesicular stomatitis virus, with minimal cytotoxicity observed in rabbit lung cells (RL-33) [1].
Targets(IC50)	Others,Antibiotic,Influenza Virus

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5357 mL	12.6784 mL	25.3569 mL
5 mM	0.5071 mL	2.5357 mL	5.0714 mL
10 mM	0.2536 mL	1.2678 mL	2.5357 mL
50 mM	0.0507 mL	0.2536 mL	0.5071 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481