

Loxosceles Spiniceps Venom

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Loxosceles Spiniceps Venom, derived from Loxosceles Spiniceps, is a complex venom composed of proteins, glycoproteins, and low molecular mass peptides. This venom contains various protein classes such as astacin-like metalloprotease, insecticidal peptides of low molecular mass, inhibitor cystine knot peptides (ICK), and phospholipase-D. The presence of the enzyme sphingomyelinase D in the venom is responsible for inducing necrosis [1].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481