


(+)-Oxypeucedanin methanolate

Chemical Properties

CAS No. :	52939-12-5
Formula:	C17H18O6
Molecular Weight:	318.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(+)-Oxypeucedanin methanolate (compound 9) is a natural inhibitor of prostaglandin E2 production [1].
Targets(IC50)	Others,Prostaglandin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1415 mL	15.7075 mL	31.4149 mL
5 mM	0.6283 mL	3.1415 mL	6.283 mL
10 mM	0.3141 mL	1.5707 mL	3.1415 mL
50 mM	0.0628 mL	0.3141 mL	0.6283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481