


Pedunculagin

Chemical Properties

CAS No. :	113866-64-1
Formula:	C ₄₁ H ₂₈ O ₂₆
Molecular Weight:	936.65
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pedunculagin, a potent 5 α -reductase type 1 inhibitor, effectively inhibits the production of nitric oxide (NO), IL-6, and IL-8 while also decreasing the protein expression of 5 α -reductase. Additionally, it exhibits anti-inflammatory activity [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0676 mL	5.3382 mL	10.6763 mL
5 mM	0.2135 mL	1.0676 mL	2.1353 mL
10 mM	0.1068 mL	0.5338 mL	1.0676 mL
50 mM	0.0214 mL	0.1068 mL	0.2135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481