


1-Stearoyl-2-arachidonoyl-sn-glycero-3-phosphoinositol

Chemical Properties

CAS No. :	1331751-28-0
Formula:	C47H83O13P
Molecular Weight:	887.13
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Stearoyl-2-arachidonoyl-sn-glycero-3-phosphoinositol is an intrinsic metabolite and a component of the cell membrane [1].
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1272 mL	5.6362 mL	11.2723 mL
5 mM	0.2254 mL	1.1272 mL	2.2545 mL
10 mM	0.1127 mL	0.5636 mL	1.1272 mL
50 mM	0.0225 mL	0.1127 mL	0.2254 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481