

GTFTSDVSKQMEEEEAVRLFIEWLKNGGPSSGAPPPS

Chemical Properties

CAS No. :

Formula: 3850.31

Molecular Weight: 3850.31

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | GTFTSDVSKQMEEEEAVRLFIEWLKNGGPSSGAPPPS is an Exendin-4 peptide derivative. |
| Targets(IC50) | Glucagon Receptor |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.2597 mL | 1.2986 mL | 2.5972 mL |
| 5 mM | 0.0519 mL | 0.2597 mL | 0.5194 mL |
| 10 mM | 0.026 mL | 0.1299 mL | 0.2597 mL |
| 50 mM | 0.0052 mL | 0.026 mL | 0.0519 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481