

## Hexokinase (ScHEX1)

## Chemical Properties

CAS No. : 9001-51-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Hexokinase (ScHEX1), a glycolytic enzyme, is inhibited by n-acetylglucosamine, which causes it to detach from the mitochondrial outer membrane. This detachment activates the NLRP3 inflammasome [1].
Targets(IC50)	Others,Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481