

Phosphoglucomutase

Chemical Properties

CAS No. : 9001-81-4

Formula:

Molecular Weight:

Storage: Store at low temperature,Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Phosphoglucomutase promotes cyst development in Toxoplasma gondii. Deficiency of Phosphoglucomutase causes metabolic disorders.
Targets(IC50)	Endogenous Metabolite

Reference

Bi-Ying Pan, et al. Role of phosphoglucomutase in regulating trehalose metabolism in Nilaparvata lugens. 3 Biotech. 2020 Feb;10(2):61.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481