

## Chymotrypsinogen

## Chemical Properties

CAS No. : 9035-75-0

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Chymotrypsinogen is an inactive precursor of Chymotrypsin . Chymotrypsin is a serine protease produced by the pancreas [1] [2] .
Targets(IC50)	Serine Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481