

Phospholipase A2

Chemical Properties

CAS No. : 9001-84-7

Formula:

Molecular Weight:

Storage: Keep away from moisture, Store at low temperature,
Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Phospholipase A2 (PLA2) is an enzyme that hydrolyzes fatty acids present in the sn-2 position of phospholipids, and is closely related to lipid metabolism and immune function.
Targets(IC50)	Phospholipase

Reference

Murakami M, Kudo I. Phospholipase A2. J Biochem. 2002 Mar;131(3):285-92.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481