

Efinopegdutide

Chemical Properties

CAS No. : 2055640-93-0

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Efinopegdutide

Biological Description

Description	Efinopegdutide (JNJ-64565111) is a potent dual agonist for the glucagon-like peptide-1 (GLP-1) and glucagon receptors (GluR), demonstrating efficacy in activating both receptors, with potential applications in the research of nonalcoholic steatohepatitis (NASH) [1].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481