


## [Asn18] Endothelin-1, human

## Chemical Properties

CAS No. :	309245-97-4
Formula:	C109H164N26O31S5
Molecular Weight:	2494.96
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	[Asn18] Endothelin-1, human, a peptide, represents a predominant isoform of Endothelin, significantly contributing to the regulation of vascular function [1].
Targets(IC50)	Endothelin Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4008 mL	2.004 mL	4.0081 mL
5 mM	0.0802 mL	0.4008 mL	0.8016 mL
10 mM	0.0401 mL	0.2004 mL	0.4008 mL
50 mM	0.008 mL	0.0401 mL	0.0802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481